ARTICULADE

SODIO IALURONATO 20 mg/ml

SODIO IALURONATO 16 mg/ml

SODIO IALURONATO 10 mg/ml

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SODIUM HYALURONATE 16 mg/ml

SODIUM HYALURONATE 10 mg/ml

VISCOELASTIC STERILE PYROGEN-FREE GEL FOR INJECTION FOR INTRARTICULAR USE

QUALITATIVE AND QUANTITATIVE COMPOSITION

SODIUM HYALURONATE 20 mg/ml

Sodium hyaluronate 25 g - Sodium chloride 8.5 g - Disodium phosphate dihydrate 0.28 g – Sodium dihydrogen phosphate dihydrate 0.04 g – Water for injections q.s. to 1000 ml.

SODIUM HYALURONATE 16 mg/ml

Sodium hyaluronate 16 g - Sodium chloride 8.5 g - Disodium phosphate dihydrate 0.28 g – Sodium dihydrogen phosphate dihydrate 0.04 g – Water for injections q.s. to 1000 ml.

SODIUM HYALURONATE 10 mg/ml

Sodium hyaluronate 10 g - Sodium chloride 8.5 g - Disodium phosphate dihydrate 0.28 g – Sodium dihydrogen phosphate dihydrate 0.04 g – Water for injections q.s. to 1000 ml.

TYPE OF DEVICE

Class III medical device. Pre-filled syringe containing a viscoelastic sterile pyrogen-free gel made of sodium hyaluronate in buffered solution. Sodium hyaluronate is obtained by fermentation and it is a high molecular weight substance (between 2.0 and 4.5 x 10⁶ Dalton).

INTENDED USE

Pain or reduced mobility due to degenerative, post-traumatic affections or changes in the joints.

ARTICULADE medical device administered intraarticularly is indicated in those patients who have had a cartilage loss and an alteration of the synovial fluid. In cases of degenerative or mechanical pathology, the administration of sodium hyaluronate temporarily restores the viscoelastic properties of synovial fluid and joint mobility with pain relief. All joints can be treated with viscosupplementation infiltrations (sodium hyaluronate), in particular the large joints that support considerable loads.

Pre-filled syringes should be used for one administration only. The injections may be repeated, but the opportunity and frequency with which the treatment can be repeated should be evaluated by the doctor, for each individual patient.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

- Do not administer intravenously
- Do not use during pregnancy and breastfeeding
- Do not use if the gel is not clear or if the primary container is damaged
- Use immediately after opening the container. The device is intended for a single and uninterrupted administration, any residue cannot be used and should be discarded due to the risk of contamination because of the loss of sterility
- Do not re-sterilize
- The device should be used only by qualified healthcare personnel, operating with aseptic technique in compliance with the relevant operating protocol. The device should be injected only once for each treatment cycle. If necessary, further injections may be given. The opportunity and frequency with

which the treatment cycle can be repeated should be evaluated by the doctor for each individual patient, considering in any case the risk/benefit ratio of the treatment.

- Do not use in patients who have had a thrombolytic or anticoagulant treatment in the previous 2 weeks
- Do not use in patients with symptoms of viral or bacterial infections or similar symptoms
- Do not use in the presence of an abundant intraarticular effusion
- After an intra-articular injection, advise the patient to avoid all intense physical activities and to resume normal activities only after a few days.

CONTRAINDICATIONS

The device has a good tolerability and acts only at the level of the area in which it is injected without exerting systemic actions.

Do not administer to patients with established individual sensitivity to sodium hyaluronate or other components present in the product.

The device should not be injected:

- If the joint is infected or severely inflamed;
- If the insertion point of the joint to be treated is infected or in the presence of a skin disease.

INTERACTIONS

In order to prevent any possible interaction, ARTICULADE should not be administered concomitantly with other products.

Sodium hyaluronate may precipitate in the presence of disinfectants containing quaternary ammonium salts or chlorhexidine.

ADVERSE REACTIONS

Localized pain, swelling, heat and redness may occur occasionally at the injection site. These symptoms may be attenuated by applying ice to the treated area and they are generally mild and transient.

More marked inflammatory reactions, sometimes with the presence of sodium pyrophosphate crystals, have occasionally been reported concomitantly with intra-articular injections of sodium hyaluronate.

As with any intra-articular treatment, septic arthritis may rarely occur when not working in asepsis.

INSTRUCTIONS FOR USE

Do not inject outside the joint cavity.

The product should be administered by qualified medical staff, adequately prepared in relation to the protocol to be applied.

Depending on the evaluation of the medical staff, needles of appropriate length and suitable caliber should be used (18-25 Gauge with Luer-Lock connector is recommended).

Screw the needle firmly to the luer-lock syringe locking collar to prevent the gel from escaping during administration. Administer sodium hyaluronate solution at room temperature and operate in compliance with aseptic conditions.

WARNINGS AND PRECAUTIONS IN RELATION TO THE SHELF LIFE OF THE DEVICE

Check the expiry date stated on the device. The expiry date refers to the product in intact packaging.

CAUTION: do not use the device after the expiry date indicated. Even within its shelf life (device not expired), the solution should not be used, if there are visible particles, precipitates or contamination.

After the first opening: after opening, the device should be used immediately for a single and uninterrupted administration. Any residues should be discarded.

SPECIAL PRECAUTIONS FOR STORAGE

The device is sterile, apyrogenic and disposable.

Store in the original package away from light at temperatures between 2° C and 25° C.

Do not freeze. Do not expose to heat sources.

KEEP OUT OF THE REACH AND SIGHT OF CHILDREN.

NATURE AND CONTENT OF THE PACK

Packaging: cardboard box showing batch number and expiry date.

Labeling: in compliance with current legislation.

Contents: Polypropylene (PP) syringe of 2 ml, 2.5 ml, 4 ml, 5 ml. Glass syringe of 2 ml, 2.5 ml, 4 ml, 5 ml.

SPECIAL PRECAUTIONS FOR DISPOSAL AND HANDLING OF THE DEVICE

The unused device and the waste derived from it should be disposed of in compliance with the local regulations in force.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP): the document is available at the following URL: <https://cloud.salfspa.it:8443/index.php/s/SSCP.Articulade>

SERIOUS INCIDENT REPORTING: the user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

MANUFACTURER: S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO Via Marconi, 2 24069 Cenate Sotto (Bergamo) – Italia

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EN – Sterile solution – The exterior part of the primary container is not sterile

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EN – Do not reuse – Use on a single patient during a single procedure

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EN – Expiry date – Date after which the medical device is not to be used

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EN – CE Mark – the product complies with applicable European regulations

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EN - The volumetric marks reported on the syringes are indicative and do not have a measurement value.