SERVATOR® B PLUS

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DATA DI REVISIONE DEL TESTO: 05/2025

SERVATOR® B PLUS

SOLUTION FOR ORGAN FLUSHING AND PRESERVATION STERILE AND PYROGEN-FREE
NOT INJECTABLE

DESCRIPTION: SERVATOR® B Plus is used for hypothermic flushing and storage of organs including kidney, liver and pancreas for transportation prior to transplantation.

COMPOSITION: 1000 ml of solution contain:

PEG 35 Kd	1.05 g/L	0.03 mmol/ L
Lactobionic acid	35.83 g/L	100 mmol/L
Sodium Hydroxide	4.8 g/L	120 mmol/L
Adenosine	1.34 g/L	5 mmol/L
Allopurinol	0.136 g/L	1 mmol/L
Potassium Phosphate monobasic	3.4 g/L	25 mmol/L
Magnesium Sulfate heptahydrate	1.23 g/L	5 mmol/L
Raffinose pentahydrate	17.83 g/L	30 mmol/L
Glutathione	0.922 g/L	3 mmol/L
Sodium hydroxide/Phosphoric acid to pH adjustment q.s. to pH 7.4		
Water for injection q.s. to 1000 mL		

PHYSICAL PROPERTIES: The solution has an osmolarity of 290 mOsm, a sodium concentration of 120 mmol/L, a potassium concentration of 25 mmol/L and a pH of 7,1 – 7,5 at room temperature.

TYPE OF DEVICE: Sterile and pyrogen-free solution of Class III for organ flushing. It is clear, colourless or slightly yellow solution.

INTENDED USE: SERVATOR® B Plus is intended for flushing and cold storage of organs including kidney, liver and pancreas at the time of their removal from the donor in preparation for storage, transportation and eventual transplantation into a recipient. It provides better preservation by reducing ischaemic-reperfusion injuries. After precooling the solution to about 2-6°C, the cold solution is used to flush the isolated organ immediately before removal from the donor and/or immediately after removal from the donor. The solution is then left in the organ vasculature during hypothermic storage and transportation.

The use of the solution at the recommended temperatures will effectively cool the organ and should reduce its metabolic requirements.

The device should be used only by medical staff properly trained according to established operating protocols.

CONTRA-INDICATIONS: Hypersensitivity to one of its components.

WARNING AND PRECAUTIONS: NOT FOR DIRECT INJECTION OR INTRAVENOUS INFUSION. USE EXCLUSIVELY FOR SOLID ORGANS FLUSHING AND COLD STORAGE.

Systemic administration may cause side effects that can endanger the patient's life.

CAUTION: Before reperfusion is established in the recipient, the donor organ must be washed using a physiological solution to remove any residual cold storage solution in order to prevent occurrence in the recipient of potentially serious cardiovascular complications. The incomplete elimination of residues of the preservation solution may result in hypersensitivity to any of its components.

ADVERSE REACTIONS: There are no known adverse reactions when used as directed.

If the cold storage solution is not completely eliminated from the organ after flushing with a physiological solution and prior to transplantation, the recipient may undergo cardiovascular complications, like bradyarrhythmia.

PREPARATION AND ADMINISTRATION: Precool the SERVATOR® B Plus to 2-6°C (36°- 43°F). Remove the overwrap before use. Do not use if the bag under pressure leaks or if visible particles, precipitate or contamination are detected upon visual inspection. Do not filter before use (the solution undergoes 0.2 micron filtration during production process). Use an administration kit to flush the organ with the SERVATOR® B Plus solution. Flushing should be continued until the organ is uniformly pale and the effluent is relatively clear.

Suggested minimum volumes:

- For in situ aortic flush
 - 2-4 L for adults
 - 50 mL/kg for infants
- For ex vivo infusion of liver
 - 1200 mL for adults
 - 50 mL/kg for infants
- For ex vivo infusion of pancreas
 - 300-500 mL for adults
 - 150-250 mL for infants

After recovery, the organ must be stored in a sterile container, covered completely by the cold solution and sealed aseptically and remain sterile.

The organ container should then be placed inside a second container and surrounded with ice for transportation. In any case ice must not come into direct contact with the organ.

FLUSHING

Prior to transplantation into the recipient, the organ must be flushed free of residual SERVATOR® B Plus solution. Flushing is to be conducted in compliance with the transplantation center's procedures. The operations must be implemented by medically qualified personnel.

WARNINGS AND PRECAUTIONS IN RELATION TO THE SHELF LIFE OF THE DEVICE: Check the expiration date on the container. Shelf life refers to the product in its intact packaging and correctly stored.

CAUTION: Do not use after the expiration date. Do not use the solution even before expiration if visible particles, precipitate or contamination are detected.

SPECIAL PRECAUTIONS FOR STORAGE: Sterile and disposable device, for one single and continuous administration. Discard any residue to avoid risk of contamination due to loss of sterility.

Store at a temperature of 2°C - 25°C (36°F - 77°F). Avoid exposure to heat. Do not freeze. Do not use if the solution is frozen. Store in its original container. Do not remove the overwrap until immediately prior to use. KEEP OUT OF THE REACH AND SIGHT OF CHILDREN.

SPECIAL PRECAUTIONS FOR DISPOSAL OF THE DEVICE: Unused residues of the solution and waste materials should be disposed of in conformity with the local regulations.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP): the document is available at the following URL: <https://cloud.salfspa.it:8443/index.php/s/SSCP.ServatorH>

SERIOUS INCIDENT REPORTING: the user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

MANUFACTURER: S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO via Marconi, 2 - 24069 Cenate Sotto (BG) Italy

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FABRICANTE: S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO via Marconi, 2 - 24069 Cenate Sotto (BG) Italia