ANTICOAGULANTE ACD-A S.A.L.F.
APYROGENIC STERILE SOLUTION BLOOD PRESERVATIVE
NON-INJECTABLE

QUALITATIVE AND QUANTITATIVE COMPOSITION

Glucose monohydrate 24.5 g - Citric acid monohydrate 8 g - Sodium citrate dihydrate 22 g - Water for injectable preparations, sufficient for 1000 ml.

PHYSICAL PROPERTIES

Limpid, colorless or light yellow solution with no visible particles.
pH of solution: 4.7 - 5.3

INTENDED USE

Medical device, sterile by steam Sterilization. Auxiliary products for blood transfusion.
Since it is not toxic, it is used to stop coagulation of the blood intended for transfusions.
The ACD anticoagulant (Acid Citrate Dextrose) is used to collect and preserve extracted blood samples and for transfusion procedures, such as plasmapheresis or Therapeutic Plasma Exchange (TPE). The citrate contained in the ACD anticoagulant takes calcium from the blood and makes it insoluble, thus preventing coagulation.

CONTRAINDICATIONS

Hypersensitivity to any of the components of the medical device - Severe renal insufficiency with oliguria, azotemia or anuria - Addison's disease – Hyperkalemic periodic paralysis – Acute dehydration - Cramps - Severe heart damage – Low sodium diet.

WARNINGS AND PRECAUTIONS

Do not inject intravenously but use only for the preparation of human blood. Use immediately after opening the container. The solution must be limpid, colorless or light yellow and with no visible particles.
It is for one single, uninterrupted use. Any residue must not be used but be disposed of.
Because of the lack of experience in pediatric apheresis procedures, it is strongly recommended to perform each treatment in co-operation with pediatricians and pediatric caregivers who should be available on demand during apheresis.
Do not use if the container is damaged or accidentally opened before use due to loss of sterility.

HOW TO USE

Use in combination with standard Blood collection set.
15 ml are used for 100 ml of collected blood. Blood collected in this way can be kept for up to 21 days.

INTERACTION WITH OTHER SUBSTANCES

There are no interactions when the device is used as indicated.
Drugs must not be added to this device.
The ACD anticoagulant can interact with the following drugs: dextran; oral anticoagulants; other drugs affecting the formation and functions of the platelets; ACE inhibitors, products containing aluminum (basic aluminum carbonate, aluminum hydroxide, aluminum phosphate, dihydroxy aluminum aminoacetate, dihydroxy aluminum sodium carbonate, magaldrate) since the simultaneous administration of a solution containing citrate can increase the absorption of aluminum, causing toxicity and leading to encephalopathy. Encephalopathy can occur particularly in elderly patients with chronic renal insufficiency.

POSSIBLE SIDE EFFECTS

Citrate toxicity may occur particularly during large volume exchanges, since citrate containing blood is returned to the donor and chelation of cations can continue in the systemic circulation. Consequently, metabolic complications ensue (such as hypocalcemia, hypomagnesemia and other electrolyte derangements) and may be accompanied by symptoms.

Typical symptoms include paresthesia (perioral or peripheral), lightheadedness, metallic taste, shivering or twitching, chest tightness, abdominal cramping, nausea and/or vomiting. Severe citrate reactions such as tetania, seizures, cardiac arrhythmias, hypernatremia, or metabolic alkalosis/acidosis are rare.

Since citrate is mainly metabolized in the liver, patients with liver disease can show a trend to increased citrate toxicity while patients with kidney dysfunction are at increased risk of metabolic alkalosis, especially during Therapeutic Plasma Exchange (TPE) with plasma as the replacement fluid.

Possible Adverse Event during Therapeutic Plasma Exchange (TPE) procedures typically consist of citrate toxicity which can be increased in patients on bisphosphonates therapy, vasovagal reactions and hypotension and hypothermia.

Paresthesia – perioral or peripheral – are among the most common citrate-induced side effects. However, severe citrate reactions such as tetania, hypernatremia, or metabolic alkalosis/acidosis are rare.

Patients with a central venous line are at risk for cardiac arrhythmia if high blood flow rates are applied as a high citrate concentration can affect the conduction system of the heart.

Other side effects that may occur are pruritus and urticaria and hypertension.

WARNINGS AND PRECAUTIONS CONCERNING THE VALIDITY OF THE DEVICE

Check the expiry date indicated on the container. The expiry date refers to the product in undamaged packaging.

WARNING: do not use the device after the indicated expiry date. Even if it is still within its expiry date, the solution must not be used if there are any visible particles, precipitates or contaminants.

SPECIAL STORAGE PRECAUTIONS

There are no special storage conditions in relation to temperature. Do not freeze.

Keep the container hermetically sealed and protected from light and heat. The device is sterile and disposable. The solution must be used for one single uninterrupted administration and any residue must be disposed of in order to prevent the risk of contamination due to loss of sterility.

KEEP OUT OF THE SIGHT AND REACH OF CHILDREN

SPECIAL PRECAUTIONS FOR THE DISPOSAL OF THE DEVICE

An unused device and the resulting waste must be disposed of in compliance with current local legislation.

Note that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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DATE OF TEXT REVISION: 11/2024

Symbols

The UNI EN ISO 15223 symbols are indicated on the packaging (Symbols to be used on the medical device labels, labeling and information that must be provided).

The symbols used and relative descriptions are listed below:

Medical Device Code

Batch number

Expiry date

The device may be broken or damaged if it is not handled with care

Protect from light

Do not reuse / Single use only. This product is exclusively disposable.

Reuse can lead to the risk of contamination and loss of sterility.

Do not re-sterilise. The materials used to produce this device could get damaged and lead to cross-contamination, deterioration, biocompatibility problems, endotoxic reactions and/or faults preventing the device from performing its proper function.

Warning: read the instructions carefully

Warning connected with a specific feature of the device

Vapor sterilisation method

Do not use if the container is damaged

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S.A.L.F. S.p.A. Laboratorio Farmacologico