Package Leaflet of the product - Notice du produit - Beipackzettel
Prospecto del producto - Prospectul produsului

(AC 113 mmol/L 2:1)
Citrate Anticoagulant Solution 113 mmol/L

Anticoagulation Solution of Trisodium Citrate (74.8 mmol/L) and Citric Acid (38.1 mmol/L)
Medical device IIb
Read carefully the following information for a safe usage of the solution.

1. Product information:

Product name:
Citrate Anticoagulation Solution 113 mmol/L (AC 113 mmol/L 2:1)
Composition:

- Each 1000 mL contains
- Natrii Citras Dihydricum Ph Eur – 22.0 g
- Acidum Citricum Monohydricum Ph Eur – 8.0 g
- Aqua ad Iniectabilia Ph Eur
- pH 4.5-6.5

Composition in mmol/L units:

Natrii Citras 74.8 mmol/L
Acidum Citricum 38.1 mmol/L

Type of device:

Citrate Anticoagulation Solution 113 mmol/L is a sterile, non-pyrogenic, clear anticoagulation solution filled in transparent printed PVC free bag and packed in polypropylene foil (overbag).
It is delivered in bags of 2000 mL.

2. Indications:

Anticoagulation solution for whole blood.
Citrate Anticoagulation Solution 113 mmol/L is used exclusively for whole blood anticoagulation within CRRT – Continuous Renal Replacement Therapy.
The solution is NOT intended for direct intravenous infusion.
The solution can be also used as a citrate lock of intravascular catheters. Citrate solution has anticoagulant properties and interferes with creation of biofilms in catheters. The citrate lock which is applied into the catheter, when no dialysis is performed, protects the catheter from thrombosis and limits creation of a catheter infection.

3. Important information:

Contraindications:

Severe liver failure and/or severe shock state with failure of citrate utilization in intermediary metabolism represent relative contraindications.

Precautions for use:

- Not to be administered intravenously.
- Not to be used unless the solution is clear.
- Check the bag and safety locks on tubing connections.
- Solution shall not be used if the container or connector safety lock has been tampered with.
- Residual fluid shall not be reused and shall be disposed of.

Interactions with medications and other types of interactions:

Citrate Anticoagulation Solution 113 mmol/L must not be added to other medications.

Warning and precautions:

- Citrate which has not been removed directly via renal hemofiltration is metabolized in the body into bicarbonate. This shall be taken into account in compounding the prescribed substitution solution.
- If Citrate Anticoagulation Solution 113 mmol/L is used during continuous blood purification procedures, it is necessary to monitor regularly all electrolytes' levels and acid base equilibrium of the patient blood.

4. Directions for use:

- Citrate Anticoagulation Solution 113 mmol/L is designed for anticoagulation of whole blood as a part of the continuous blood purification method, CRRT. The solution can be used also as a citrate lock for intravascular catheters. The citrate lock which is applied into the catheter, when no dialysis is performed, protects the catheter from thrombosis and limits creation of a catheter infection.
- Citrate Anticoagulation is considered to be a suitable method of anticoagulation mainly to use in patients with risk of hemorrhage and is also recommended to use in patients with critical conditions hospitalized in intensive care units.
- Citrate Anticoagulation Solution 113 mmol/L is infused in pre-dilution. The ratio of citrate to blood flow rate should be set to provide a final citrate concentration in blood 3-5 mmol/L. This concentration decreases ionized calcium levels inside the filter below 0.4 mmol/L and provides sufficient anticoagulation. Around 35-50% of the citrate is, depending on the configuration of CRRT, administered in the patient's systemic circulation. The concentration level of ionized calcium shall be monitored to keep the level in physiological limits. If necessary, a calcium replacement solution is administered by infusion, right before the connection to the venous catheter. A suitable starting dose is 1.7 mmol of calcium per 1 liter of the applied dialysis solution and is further adjusted to the level of ionized calcium in arterial blood which should be higher than 0.80 mmol/L. The shift to acidosis and compensation of the possible alkalizing effect of the trisodium citrate present should be taken into account. Vice versa, increase of the blood flow and related doses of the citrate can result in an increase of the citrate concentration and a shift to alkalosis. Modality, next to the regular monitoring of ionized calcium in blood and behind the filter, requires administration of a special dialysis hyponatremic solution with a lowered level of alkali to avoid a metabolic alkalosis and hypernatremia.
- Remove the bag from the cover immediately before use.

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- Check the composition, production batch number and expiry date.
- Check the solution for clarity.
- Check the bag and safety locks on tubing connectors if tamper-proof. The solution shall not be used if container or outlet safety lock has been tampered with.
- Carry out the blood purification method according to detailed operating instructions of the manufacturer of dialysis monitors. The amount of Citrate Anticoagulation Solution 113 mmol/L administered via pre-dilution shall be included in the fluid balance and the ultrafiltration calculation.
- When unpacking the bag from the secondary packaging, do not grasp tubing connectors.
- Not for re-use: Citrate Anticoagulation Solution 113 mmol/L is for single use only and is not intended to be reprocessed or reused. Reuse can lead to the risk of contamination and loss of sterility.
- Discard unused portion.

5. Adverse reactions:

Main adverse reactions are: hypocalcemia connected with tetany, paresthesia and heart function disorders (QT interval, arrhythmia, hypocontractility of myocardium). Routine substitution of calcium chloratum during the dialysis should prevent mentioned effects. It is necessary to monitor the total to ionized calcium index and assume accumulation of the citrate in organism when the index is higher than 2.5.
Other adverse reactions are: metabolic alkalosis, hypernatremia, hypomagnesemia.

6. Expiry date and storage conditions:

Expiry date of Citrate Anticoagulation Solution 113 mmol/L is imprinted on the container's rear side. The product shall not be used after this date.
Anticoagulation Citrate Solution 113 mmol/L shall be protected from light and stored between 5-25°C. Do not freeze.

Packing: 5 x 2000 mL PVC free bags.

7. Information on packagings' disposal:

The product is recommended to be considered potentially hazardous waste.
Note that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Symbols

The symbols stated on the container and relative descriptions are listed below:

	Medical Device Code
	LOT Number
	Expiry Date
	Protect from direct sun light
	Do not reuse / Single use only. This product is exclusively disposable. Reuse can lead to the risk of contamination and loss of sterility.
	Do not re-sterilise. The materials used to produce this device could get damaged and lead to cross-contamination, deterioration, biocompatibility problems, endotoxic reactions and/or faults preventing the device from performing its proper function.
	Pay attention to the instructions for use!
	Indicate that the product is non-pyrogenic
	Sterilization by steam method
	Single sterile barrier system with protective packaging outside
	Do not use if the immediate container is damaged
	Temperature limit
	Manufacturer

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S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO via Marconi, 2 – 24069 Cenate Sotto (BG) – Italia