OCULADE

SODIO IALURONATO 25 mg/ml

SODIO IALURONATO 15 mg/ml

SODIO IALURONATO 14 mg/ml

SODIO IALURONATO 10 mg/ml

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VISCOELASTIC STERILE PYROGEN-FREE GEL FOR INJECTION FOR OPHTHALMIC USE

QUALITATIVE AND QUANTITATIVE COMPOSITION

SODIUM HYALURONATE 25 mg/ml

Sodium hyaluronate 25 g - Sodium chloride 8.5 g – Disodium phosphate dihydrate 0.28 g – Sodium dihydrogen phosphate dihydrate 0.04 g – Water for injections q.s. to 1000 ml.

SODIUM HYALURONATE 15 mg/ml

Sodium hyaluronate 15 g - Sodium chloride 8.5 g – Disodium phosphate dihydrate 0.28 g – Sodium dihydrogen phosphate dihydrate 0.04 g – Water for injections q.s. to 1000 ml.

SODIUM HYALURONATE 14 mg/ml

Sodium hyaluronate 14 g - Sodium chloride 8.5 g – Disodium phosphate dihydrate 0.28 g – Sodium dihydrogen phosphate dihydrate 0.04 g – Water for injections q.s. to 1000 ml.

SODIUM HYALURONATE 10 mg/ml

Sodium hyaluronate 10 g - Sodium chloride 8.5 g – Disodium phosphate dihydrate 0.28 g – Sodium dihydrogen phosphate dihydrate 0.04 g – Water for injections q.s. to 1000 ml.

TYPE OF DEVICE

Class III medical device. Pre-filled syringe containing a viscoelastic sterile pyrogen-free gel made of sodium hyaluronate in buffered solution. Sodium hyaluronate is obtained by fermentation and it is a high molecular weight substance (between 2.0 and 4.5 x 10⁶ Dalton).

The syringe is equipped with a protective cover to guarantee the sterility of its external part.

INTENDED USE

OCULADE medical device is indicated in ophthalmic surgery, both to maintain an adequate depth of the endocular spaces of the anterior and posterior segments of the eye and to preserve the integrity of the tissues from mechanical damage during surgery.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

Do not administer intravenously.

Do not use if the protective cover is damaged, if the primary container is damaged or if the gel is not clear.

Use immediately after opening the container. The device is intended for a single and uninterrupted administration, any residue cannot be used and should be discarded due to the risk of contamination because of the loss of sterility. Do not re-sterilize.

The device should be used only by qualified healthcare personnel operating with aseptic technique in compliance with the relevant operating protocol.

CONTRAINDICATIONS

The device has a good tolerability and acts only at the level of the area in which it is injected without exerting systemic actions.

Do not administer to patients with established individual sensitivity to sodium hyaluronate or other components of the product.

The device should not be used in the presence of inflammatory processes.

INTERACTIONS

Sodium hyaluronate may precipitate in the presence of disinfectants containing quaternary ammonium salts or chlorhexidine.

ADVERSE REACTIONS AND SIDE EFFECTS

The solution of sodium hyaluronate, if not correctly removed at the end of the surgery, may cause an increase in intraocular pressure and/or a transient corneal oedema in the postoperative phase.

It is therefore necessary to follow the following precautions:

- In order to significantly reduce the increase in intra-ocular pressure, it is essential to remove excess device by aspiration and/or irrigation with sterile saline solution before closing the incision. In any case, at the end of the surgical intervention the patient's intra-ocular pressure should be carefully monitored, in order to intervene with appropriate therapy in case of need
- Pre-existing glaucoma or high preoperative intraocular pressure could lead to excessive intraocular pressure immediately after the operation. For this reason, in these cases it is of particular importance to monitor the intraocular pressure immediately after the surgical procedure.

Following surgery, rare cases of inflammatory reaction (including TASS: toxic anterior segment syndrome), corneal oedema, secondary glaucoma, IOP (intraocular pressure) spikes associated with increased risk of vision loss, corneal de-compensation, Capsular Block Syndrome, postoperative uveitis, hypopyon, IOL (intraocular lens) surface crystallization and posterior capsule rupture have been reported.

These problems have not been demonstrated as attributable to the presence of sodium hyaluronate.

INSTRUCTIONS FOR USE

The product should be administered by qualified medical staff, adequately prepared in relation to the protocol to be applied, while operating in compliance with the conditions of asepsis.

Bring the device to room temperature before use.

Depending on the doctor's evaluation, sterile and disposable cannulas of appropriate length, suitable caliber and equipped with Luer Lock connector should be used (depending on the concentration of the product, a 25 or 27 Gauge needle is recommended).

The amount of device administered by injection in the anterior chamber at different stages of surgery is related to the type of surgery. Before the incision is closed, the excess device should be removed from the anterior or posterior chamber by aspiration and/or irrigation.

WARNINGS AND PRECAUTIONS IN RELATION TO THE SHELF LIFE OF THE DEVICE

Check the expiry date stated on the primary container of the device. The expiry date refers to the product in intact packaging.

CAUTION: do not use the device after the expiry date indicated. Even within its shelf life (device not expired), the solution should not be used, if there are visible particles, precipitates or contamination.

After the first opening the device should be used immediately for a single and uninterrupted administration. Any residue should be discarded.

SPECIAL PRECAUTIONS FOR STORAGE

The device is sterile and disposable.

Store in the original package away from light at a temperature between 2° C and 25° C.

Do not freeze. Do not expose to heat sources.

Keep out of the reach and sight of children.

NATURE AND CONTENT OF THE PACK

Packaging: cardboard box showing batch number and expiry date.

Labeling: in compliance with current legislation.

Contents: polypropylene or glass syringe of 1 ml or 0.5 ml.

SPECIAL PRECAUTIONS FOR DISPOSAL AND HANDLING OF THE DEVICE

The unused device and the waste derived from it should be disposed of in compliance with the local regulations in force.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP): the document is available at the following URL:
<https://cloud.salfspa.it:8443/index.php/s/SSCP.Oculade>

SERIOUS INCIDENT REPORTING: the user and/or patient are recommended to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

MANUFACTURER: S.A.L.F. S.p.A. LABORATORIO FARMACOLOGICO Via Marconi, 2 24069 Cenate Sotto (Bergamo) – Italia

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EN – Sterile solution – The exterior part of the primary container is not sterile

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EN – Do not reuse – Use on a single patient during a single procedure

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EN – Expiry date – Date after which the medical device is not to be used

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EN – CE Mark – the product complies with applicable European regulations

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EN – Single sterile barrier system with protective packaging outside

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EN - The volumetric marks reported on the syringes are indicative and do not have a measurement value.