EN
How to use

Sodium citrate 3.8%
sterile non-pyrogenic anticoagulant solution

COMPOSITION
Sodium citrate dihydrate 38g, water for injection sufficient quantity for 1,000 ml.

CLASSIFICATION
Sterile hypertonic solution, non-pyrogenic, Class IIb.
The device is sterile for terminal sterilization with saturated steam.

USE
The solution (all or part) is used for the preparation of autologous Plasma Rich in Growth Factors and Platelet Rich Plasma (PRP) fraction as a whole blood anticoagulant before and during the separation and concentration of the blood component fractions by gradient density. Platelets contain growth factors within their "alpha" granules, which are responsible for physiological tissue repair. This therapy is widely described in particular for promoting the repair of soft tissues and bones, especially in difficult treatment cases, where the regeneration process is slow and complicated.

SHELF LIFE
5 years from date of manufacture, where the packaging is intact.

STORAGE CONDITIONS
There are no special storage conditions.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

- Do not administer intravenously.
- Before administration, make sure that the solution is clear and free from visible particles.
- Do not use if the container is damaged or accidentally opened before use.
- Use immediately after opening the ampoule. Single-use solution.
- Do not use any residual solution in any other treatment. Discard the open ampoule and only use a new sealed ampoule.

CONTRAINDICATIONS
No adverse reactions or contraindications are known, provided that the indications and directions for use are followed properly.

INSTRUCTIONS FOR USE
The device can only be used by properly trained medical personnel according to established operating protocols.

Follow the protocol for use, as indicated: prior to venipuncture and removal of peripheral venous blood, preload the suction device (syringes or specific devices) with a volume of 3.8% sodium citrate equal to 10% of total blood to be removed (1 volume of sodium citrate to 9 volumes of blood – 1:9).

INTERACTIONS: In the absence of incompatibility studies, Sodium Citrate 3.8% must not be mixed with other products.

CONTRAINDICATIONS: Do not directly infuse the product as it is into the vein.

UNDESIRABLE EFFECTS: No adverse events are known if Sodium Citrate 3.8% Solution is used correctly following the specific Apheresis protocol.

OVERDOSE: Sodium Citrate 3.8% solution is used in quantities specified by the relevant apheresis protocol and should not be administered by direct intravenous infusion.

PACKING

- Box with 5 glass ampoules containing 10 ml of solution.
- Box with 5 glass ampoules containing 5 ml of solution.
- Box with 5 glass ampoules containing 4 ml of solution.
- Box with 5 glass ampoules containing 2 ml of solution.

Note that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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